

deconex® CIP seven

Neutral special cleaning concentrate with a broad material compatibility for the pharmaceutical industry

Liquid, free of surfactants. For spray cleaning processes in CIP-systems, in washers and for manual cleaning.



Application

deconex® CIP seven is used for an effective cleaning of different materials in the pharmaceutical industry. Validated cleaning processes based on deconex® CIP seven are suitable for the

- pharmaceutical industry
- biotechnological production

A residue-free cleaning of surfaces in contact with pharmaceutical products is an important condition in critical areas such as the manufacturing of drugs. The application of deconex® CIP seven is recommended for spray cleaning in CIP-systems or in washers.

deconex® CIP seven is also efficient for flooding of production equipment, mixers or containers as well as for the manual cleaning of spare parts.

deconex® CIP seven is an excellent detergent for the efficient cleaning of peripheral parts in the processing of solid dosage forms.

Properties

deconex® CIP seven is

- pH-neutral in the ready-to-use solution
- compatible with aluminium
- suitable for anodized aluminium
- free of silicates and surfactants
- also working without neutralization step
- easily rinsed-off
- absolutely non-foaming
- suitable for validated cleaning processes

Save time by using deconex® CIP seven!

The mixed loading and reprocessing in a washer of objects made from different materials such as aluminium, plastics, stainless steel and glass is absolutely possible.

Ingredients

Alkalis, complexing agents

Dosage

The dosage depends on the type of residues, the substrates, the process and the water hardness.

Usually, the optimum concentration in use is 0.5-2.0 % (v/v).

Information on use

The addition of deconex® CIP surf for the elimination of particularly stubborn residues shows excellent results.

Neutralization

A neutralization step with e.g. deconex® CIP fresh may shorten the rinsing steps and rinsing time considerably. If cleaning is performed in tap water the neutralization step eliminates possible residues from mineral and metal ions.

Residue analysis/cleaning validation

In order to confirm a residue-free cleaning suitable and validated analytical methods are supplied. We are pleased to offer advice on how to perform them in practice.

deconex® CIP seven

Material compatibility

Suitable for:
Aluminium, anodized aluminium, stainless steel, plastics, glass

For materials not mentioned please make your own specific compatibility tests or consult Borer Chemie AG.

Chemical/physical data

pH	1% solution in demineralized water	approx. 8.8
Density	concentrate	1.2 g/ml
Appearance	concentrate	transparent, light yellowish

Availability

Please ask your local representative about current container sizes.

Containers, screw caps and labels are made of recyclable polyethylene.

Additional information

For information concerning safety at work, storage and waste disposal/effluent, please consult the corresponding safety data sheet.

Take advantage of our vast know-how! Please, contact us for further information regarding your specific application.

deconex® CIP for the highest standards of quality

deconex® CIP products have been developed specially for validated cleaning.

Our products and services allow individual, optimally adapted and efficient cleaning processes to be used. You can obtain further information on this directly from Borer Chemie AG.

Manufacturer:

Borer Chemie AG

Gewerbestrasse 13, 4528 Zuchwil / Switzerland
Tel +41 32 686 56 00 Fax +41 32 686 56 90
office@borer.ch, www.borer.ch

All information provided is based on our current knowledge and it does not constitute a legally binding assurance of specific product properties.



advanced cleaning solutions