

deconex® 32 EMULGATOR

Liquid additif for the emulsification of fatty substances

For use in laboratory washing machines



Application

deconex 32 EMULGATOR is used as an additive to acid or alkaline cleaning agents and to emulsify organic residues like greases and oils. It is applied in the pre- and main-rinsing cycle.

Properties

deconex 32 EMULGATOR is a liquid, neutral, highly active and weakly foaming emulsifier. Moreover, deconex 32 EMUL-

GATOR is like a wetting agent with a foam reducing effect. It leaves a nicely shiny effect on surfaces and it is at 97% biodegradable.

Ingredients

Surfactant compounds

Dosage

Rinsing solution	0.5-2 ml/l
------------------	------------

Information:

If possible, add dosage only at a temperature of 40 °C. Please observe the recommendations of the washing machine manufacturer.

Chemical/physical data

pH	concentrate	approx. 6.5 +/-0.5
Density	concentrate	1.01 g/ml
Appearance	concentrate	clear, colourless liquid

Availability

Please ask your local representative about current container sizes.

Containers, screw caps and labels are made of recyclable polyethylene.

Additional information

For information concerning safety at work, storage and waste disposal/effluent, please consult the corresponding safety data sheet.

Take advantage of our vast know-how! Please, contact us for further information regarding your specific application.

deconex® 32 EMULGATOR

Manufacturer:

Borer Chemie AG

Gewerbestrasse 13, 4528 Zuchwil / Switzerland

Tel +41 32 686 56 00 Fax +41 32 686 56 90

office@borer.ch, www.borer.ch

All information provided is based on our current knowledge and it does not constitute a legally binding assurance of specific product properties.